

The iBioSim User's Manual

Chris J. Myers, Nathan Barker, Hiroyuki Kuwahara, Curtis Madsen, Nam Nguyen

Created: August 6th, 2008
Last Revised: August 6th, 2008

Contents

| | | |
|----------|--|-----------|
| 1 | Introduction | 3 |
| 2 | Installation Instructions | 3 |
| 2.1 | Installation on Windows | 3 |
| 2.2 | Installation on Linux | 3 |
| 2.3 | Installation on MacOS | 3 |
| 3 | Project Management | 4 |
| 3.1 | Creating and Opening Projects | 4 |
| 3.2 | Creating Models and Graphs | 5 |
| 3.3 | Importing Models | 5 |
| 3.4 | Creating and Opening Views | 6 |
| 3.5 | Viewing and Modifying Project Objects | 6 |
| 4 | SBML Editor | 7 |
| 4.1 | SBML Math Formulas | 8 |
| 4.2 | Main Elements | 9 |
| 4.2.1 | Compartments | 9 |
| 4.2.2 | Species | 10 |
| 4.2.3 | Reactions | 11 |
| 4.2.4 | Global Parameters | 12 |
| 4.3 | Definitions/Types | 13 |
| 4.3.1 | Function Definitions | 13 |
| 4.3.2 | Unit Definitions | 14 |
| 4.3.3 | Compartment Types | 14 |
| 4.3.4 | Species Types | 15 |
| 4.4 | Initial Assignments/Rules/Constraints/Events | 15 |
| 4.4.1 | Initial Assignments | 15 |
| 4.4.2 | Rules | 16 |
| 4.4.3 | Constraints | 16 |
| 4.4.4 | Events | 17 |
| 5 | GCM Editor | 18 |
| 5.1 | Promoters | 19 |
| 5.2 | GCM Species | 19 |
| 5.3 | Influences | 20 |
| 5.4 | GCM Parameters | 21 |
| 6 | Analysis View | 22 |
| 6.1 | Simulation Options | 22 |
| 6.2 | Abstraction Options | 23 |
| 6.3 | Parameter Editor | 24 |
| 7 | Learn View | 25 |
| 7.1 | Data Manager | 25 |
| 7.2 | Learn Tool | 26 |

| | | |
|-----------|--|-----------|
| 8 | TSD Graph Editor | 28 |
| 9 | Probability Graph Editor | 30 |
| 10 | Genetic Circuit Model Format | 32 |
| 11 | Time Series Data Format | 33 |
| 12 | Tutorial | 34 |
| 13 | Reporting Bugs and Feature Requests | 38 |
| 14 | Credits | 38 |

1 Introduction

iBioSim has been developed for the analysis of biochemical reaction network models. While the primary target of iBioSim is models of genetic circuits, models representing metabolic networks, cell-signaling pathways, and other biological and chemical systems can also be analyzed. iBioSim includes the following tools:

- GCM Editor - a tool to create a model using the *Genetic Circuit Model* (GCM) format. The GCM format is described in Nguyen's MS Thesis (UofUtah 2008).
- SBML Editor - a tool to create a model using the *Systems Biology Markup Language* (SBML).
- reb2sac - an abstraction-based ODE, Monte Carlo, and Markov analysis tool. This tool is described in Kuwahara's PhD Dissertation (UofUtah 2007).
- GeneNet - a tool to learn a GCM from *time series data* (TSD). This tool is described in Barker's PhD Dissertation (UofUtah 2007).
- TSD Graph Editor- a tool to visualize TSD files.
- Probability Graph Editor - a tool to visualize probability data.

2 Installation Instructions

There are versions of iBioSim available for Windows, Linux, and MacOS. You can download the appropriate installation file from:

<http://www.async.ece.utah.edu/iBioSim> To install follow the instructions for your operating system below.

DESCRIBE Java, Graphviz requirements

2.1 Installation on Windows

Download and execute iBioSim-`<version>`-Setup.exe. MORE HERE

EXPLAIN UNINSTALL

2.2 Installation on Linux

Download iBioSim-`<version>`-Linux-x86-Install. Change permissions on this file to allow it to be executed:

`chmod u+x iBioSim-<version>-Linux-x86-Install`. Executing this file will now start the install shield.

EXPLAIN UNINSTALL

2.3 Installation on MacOS

EXPLAIN IT HERE

EXPLAIN UNINSTALL


3 Project Management

A project is a collection of models, analysis views, learn views, and graphs. As shown below, iBioSim displays all project files on the left, the open models, views, and graphs are on the right, and a log of all external commands are displayed on the bottom. The menu bar is located on the top of the window in the Windows and Linux versions. It is located on the top of the screen in the MacOS version.




3.1 Creating and Opening Projects

To create a new project, select New → Project from the File menu as shown below. You will then be prompted to browse to a desired location and to give a name to the project directory. After you do this, click the new button and a new project directory will be created. To open a project, select Open → Project from the File menu. You will then be prompted to browse to a project directory to open, and clicking open will open the project. You may also open a project by selecting one of your five most recently opened projects by selecting the project name shown in the File drop down menu shown below.




3.2 Creating Models and Graphs

After you have created or opened a project, you can create a new model or graph to add to the project. To create a new Genetic Circuit Model (see Section 10), select New → Genetic Circuit Model from the File menu as shown below. You will then be prompted to give a model id. At this point, a GCM editor (see Section 5) will open in a new tab. To create a new SBML model, select New → SBML Model from the File menu. You will then be prompted to give a model id. At this point, an SBML editor (see Section 4) will open in a new tab. To create a new TSD graph, select New → TSD Graph from the File menu. You will then be prompted to give a name to the TSD graph. At this point, a TSD graph editor (see Section 8) will open in a new tab. To create a new probability graph, select New → Probability Graph from the File menu. You will then be prompted to give a name to the probability graph. At this point, a probability graph editor (see Section 9) will open in a new tab. Once a model or graph is created, it can be opened again later by right clicking on the object in the project window and selecting “Edit”, or alternatively double-clicking on the object.




3.3 Importing Models

You can import into the current working project GCMs or SBML Models created by other programs or stored in other projects. To import a GCM, select Import → Genetic Circuit Model from the File menu as shown below. You will then be able to browse to find a model to import. After selecting the desired model, click the import button to bring the GCM into the project. Before bringing the model into the project, it will be checked to see if it is a valid GCM file. To import an SBML model, the procedure is the same except use the Import → SBML Model option. Before bringing the model into the project, it will be checked to see if it is a valid SBML file. The model will also be checked for consistency, and any errors or warnings will be reported. These should be corrected before analysis of the model is performed.



3.4 Creating and Opening Views

To perform analysis or learning, right click on a model and select “Create Analysis View” (see Section 6) to perform analysis or “Create Learn View” (see Section 7) to perform learning. You will then be prompted to give a name to your analysis or learn view. After a name is entered, a tab with the newly created view will open. Once a view is created, it can be opened again later by right clicking on an analysis directory and selecting “Open Analysis/Learn View” or alternatively double-clicking on the view.




When you create an analysis view from a GCM, an SBML model is automatically created for simulation and analysis. Within the analysis view, you can edit the initial concentrations and parameters. However, if you wish to be able to edit the structure, you should first create an SBML model using Create SBML Model option in the right click menu or the Save as SBML button in the GCM Editor (see Section 5). You can then open and edit this model using an SBML editor (see Section 4) and create an analysis view from this edited model.

3.5 Viewing and Modifying Project Objects


All project objects can be modified by highlighting the object and using a right mouse click to open a menu of options as shown below. Using this menu, every type of object can be copied, renamed, or deleted. For a GCM, the “View/Edit” option opens the model in a GCM editor (see Section 5). For an SBML model, the “View/Edit” option opens the model in a SBML editor (see Section 4). For a TSD graph, the “View/Edit” option opens the TSD graph in a TSD graph editor (see Section 8). For a probability graph, the “View/Edit” option opens the probability graph in a probability graph editor (see Section 9). For an analysis view, the “Open Analysis View” option opens the analysis view (see Section 6). For a learn view, the “Open Learn View” option opens the learn view (see Section 7).

A GCM can also be viewed using GraphViz’s `dot` program by right clicking on the model you want to view and selecting the “View Genetic Circuit” option. There are two additional ways to view an SBML model. You can either select the “View Network” option or the “View in Browser” option. The “View Network” option converts the model to a GraphViz file and then will open that file with GraphViz’s `dot` program. The “View in Browser” option converts the model to an xhtml file and opens that file with your default xhtml browser.



4 SBML Editor

The SBML editor as shown below allows the user to create or modify an SBML model of a biochemical reaction network. An SBML model includes compartments (see Section 4.2.1), species (see Section 4.2.2), reactions (see Section 4.2.3), parameters (see Section 4.2.4), function definitions (see Section 4.3.1), unit definitions (see Section 4.3.2), compartment types (see Section 4.3.3), species types (see Section 4.3.4), initial assignments (see Section 4.4.1), rules (see Section 4.4.2), constraints (see Section 4.4.3), and events (see Section 4.4.4). Each of these items can be added, removed, or edited. To add a new item, click on the appropriate add button. You will then be prompted to provide a unique id and some properties for this new item (as described below). After you have filled out all of the required fields, click add and the new item will be added to the SBML model. To remove an item from the model, select that item and click the remove button. The item will then be removed from the model. However, if you try to remove an item that is being used (for example, a species that is used in a reaction), you will first have to remove its use. To edit an existing item, select that item from the list and click the edit button. An editing window will open and you will be able to change the properties of that item. When you are done editing this item, click save to save the changes to the item. After the model is complete, press the Save SBML button to store the model. The Save and Check SBML button also saves the model, but in this case it also checks the models consistency. Finally, the Save As button can also be used to store the model, but in this case, a new model ID will be requested and the model will be saved using that ID.



4.1 SBML Math Formulas

Math formulas appear in many SBML constructs. These formulas are expressed as text strings using a simple C-like syntax. SBML math formulas can include:

- Variables (compartment, species, parameter IDs, and reaction IDs)
- Real Numbers
- Built-in constants: exponentiale, pi, true, and false.
- Special variable time or t which returns the current simulation time.
- Mathematical operators including add (+), subtract (-), multiply (*), divide (/), and power (x^y) which is equivalent to pow(x,y).
- A function defined in the list of function definitions.
- Logical functions: and, or, xor, not.
- Relational functions: eq, neq, geq, gt, leq, and lt.
- Unary functions: abs, ceiling, exp, factorial, floor, ln, log, sqr, and sqrt.
- Trigonometric functions: cos, cosh, sin, sinh, tan, tanh, cot, coth, csc, csch, sec, sech, arccos, arccosh, arcsin, arcsinh, arctan, arctanh, arccot, arccoth, arccsc, arccsch, arcsec, and arcsech.
- The delay(expr1,expr2) function which returns the value of expr1 at a time expr2 time units earlier (not currently supported by analysis).
- The piecewise(value1, case1, value2, case2, ..., otherwise) function returns value1 if case1 is true, value2 if case2 is true, etc. If no cases are true, it returns otherwise value.

iBioSim’s simulators also support several random functions which are added by default to any SBML file created with its SBML editor. The following random functions, therefore, can also be used in SBML math formula:

- Continuous random functions: `uniform(a,b)`, `normal(m,s)`, `exponential(mu)`, `gamma(a,b)`, `log-normal(z,s)`, `chisq(nu)`, `laplace(a)`, `cauchy(a)`, and `rayleigh(s)`.
- Discrete random functions: `poisson(mu)`, `binomial(p,n)`, and `bernoulli(p)`.

4.2 Main Elements


This main elements tab shown below is used to specify compartments (see Section 4.2.1), species (see Section 4.2.2), reactions (see Section 4.2.3), and parameters (see Section 4.2.4). This tab also includes the Model ID which is fixed to be the same as the filename as well as the Model Name which can be used to provide an arbitrary string description of the model.

4.2.1 Compartments

Compartments are used to specify locations where species are found. Every model must include at least one compartment. A new model includes a compartment named “default” that cannot be removed unless a new compartment is provided. A compartment to which species have been assigned also cannot be removed. As shown below, a compartment has the following fields:

- ID: a unique ID composed of only alphanumeric characters and underscores.
- Name: an arbitrary string description (optional).
- Type: selected from the list of compartment types (default=none).
- Dimensions: number of spatial dimensions (default=3).
- Outside: the compartment that is outside this compartment (default=none).
- Constant: Boolean indicating if the size is constant (default=true).


- Size: initial size of the compartment (default=1.0).
- Units: the units for the size (default=none).



4.2.2 Species

Species are the molecules that appear as reactants, products, or modifiers in the reactions in the biochemical reaction network. As shown below, a species has the following fields:

- ID: a unique ID composed of only alphanumeric characters and underscores.
- Name: an arbitrary string description (optional).
- Type: selected from the list of species types (default=none).
- Compartment: location of the species (default=default).
- Boundary Condition: Boolean indicating if the species amount/concentration cannot be changed by reactions (default=false).
- Constant: Boolean indicating if the species amount/concentration is constant (default=false).
- Initial Amount/Concentration: initial value of the amount or concentration of the species whether it is an amount or concentration can also be selected (default=amount/0.0).
- Units: the units for the amount/concentration (default=none).




4.2.3 Reactions

Reactions are used to create or destroy molecular species in a biochemical reaction network. As shown below, a reaction is composed of the following:

- ID: a unique ID composed only of alphanumeric characters and underscores.
- Name: an arbitrary string description (optional).
- Reversible: a Boolean indicating if the reaction is reversible (default=false).
- Fast: a Boolean indicating if the reaction is fast (default=false).
- List of Reactants: species that are consumed by this reaction.
- List of Products: species that are produced by this reaction.
- List of Modifiers: species that are neither produced or consumed by this reaction.
- List of Local Parameters: symbolic values that can be used in the kinetic law or stoichiometry math formulas for this reaction.
- Kinetic law: a SBML math formula (see Section 4.1) describing the rate or probability for this reaction.

When adding a reactant or product, the user must specify a species ID and the stoichiometry (i.e., the number of molecules produced or consumed by the reaction). The stoichiometry can also be expressed as a stoichiometry math formula. Each parameter is composed of an ID, Name, Value, and Units. The list of parameters begins with a default forward reaction rate (kf) and reverse reaction rate (kr). These names and their values should likely be edited. The kinetic law can either be automatically generated using the Use Mass Action button or manually entered. The “Use Mass Action” button creates a rate law using the law of mass action. It assumes that the first parameter


in the list is the forward reaction rate and the second parameter in the list is the reverse reaction rate. The “Clear” button clears the kinetic law editor. The stoichiometry math and kinetic law formulas can only include those species that appear as reactants, products, or modifiers for the reaction.



4.2.4 Global Parameters


Global parameters are variables that can be used in SBML math formulas (see Section 4.1). As shown below, a parameter includes the following:

- ID: a unique ID composed only of alphanumeric characters and underscores.
- Name: an arbitrary string description (optional).
- Value: initial value for the parameter.
- Units: the units for the parameter value (default=none).
- Constant: Boolean indicating if the parameter value is constant (default=true).




4.3 Definitions/Types

The Definitions/Types tab shown below allows users to provide function definitions (see Section 4.3.1), unit definitions (see Section 4.3.2), compartment types (see Section 4.3.3), and species types (see Section 4.3.4).



4.3.1 Function Definitions

Function definitions are used to create user defined functions that can then be used in SBML math formulas (see Section 4.1). As shown below, function definitions include an ID, an optional name field, a comma-separated list of arguments, and its definition. The definition is a SBML math formula though it is restricted to only use variable names which are arguments to the function. While functions can call other functions, they cannot be recursive (i.e., call themselves) either directly or indirectly (i.e., through a cycle of function calls). As mentioned earlier, several random functions supported by iBioSim's simulators are added automatically.




4.3.2 Unit Definitions

Unit definitions are used to construct user-defined units which are derived from the set of base units. As shown below, a unit definition includes an ID, an optional name, and a list of units that define it. There are buttons to add, remove, and edit elements in the list of units. Each unit is composed of a kind, exponent, scale, and multiplier. The kind is selected from the list of base units in the table below:

| | | | | | |
|---------------|-------|----------|--------|-----------|-------|
| ampere | gram | katal | metre | second | watt |
| bacquerel | gray | kelvin | mole | siemens | weber |
| candela | henry | kilogram | newton | sievert | |
| coulomb | hertz | litre | ohm | steradian | |
| dimensionless | item | lumen | pascal | tesla | |
| farad | joule | lux | radian | volt | |


The exponent and scale are integers, and the multiplier is a real number that specify the relationship between the derived unit and the base unit using the relation below:

$$\text{unit} = (\text{multiplier} * 10^{\text{scale}} * \text{baseUnit})^{\text{exponent}}$$




4.3.3 Compartment Types

Compartment types are used to relate multiple compartments. As shown below, a compartment type includes an ID and an optional name field.




4.3.4 Species Types

Species types are used to relate multiple species. As shown below, a species type includes an ID and an optional name field.




4.4 Initial Assignments/Rules/Constraints/Events

This tab as shown below allows users to provide initial assignments (see Section 4.4.1), rules (see Section 4.4.2), constraints (see Section 4.4.3, and events (see Section 4.4.4).



4.4.1 Initial Assignments

Initial assignments as shown below are used to provide a SBML math formula (see Section 4.1) that is evaluated at time 0 to determine the initial value of a compartment size, species amount or concentration, or parameter. The value of this formula takes precedence over the initial value specified in the object.




4.4.2 Rules

There are three types of rules: algebraic, assignment, and rate rules which are in the following form:

| | | |
|------------|------------------------------------|------------------------|
| Algebraic | left-hand side is zero | $0 = f(W)$ |
| Assignment | left-hand side is a scalar | $x = f(W)$ |
| Rate | left-hand side is a rate-of-change | $\frac{dx}{dt} = f(W)$ |

Algebraic rules specify relationships which must be maintained (not currently supported by analysis). Assignment rules specify the value of a compartment size, species amount or concentration, or parameter in terms of a SBML math formula (see Section 4.1). A variable cannot be determined by both an assignment rule and initial assignment. Rate rules specify the rate of change of a compartment size, species amount or concentration, or parameter in terms of a SBML math formula (see Section 4.1). A variable cannot be determined by both an assignment rule and a rate rule. A species that is reactant or product of any reaction cannot be updated by either an assignment rule or rate rule.


When adding a rule, the user first selects the type of rule as shown below. This will automatically restrict the set of variables available for the left-hand side to those that are valid. The user should then select a variable, and enter a SBML math formula (see Section 4.1) for the rule. When editing a rule, the user cannot modify the rule type.



4.4.3 Constraints


Constraints are used to specify properties that should cause simulation to terminate. Our analysis method can provide histograms that show the proportion of simulations that are terminated due to each possible constraint. As shown below, each constraint is composed of an ID which is used to

identify it in these histograms, a constraint given as an SBML math formula (see Section 4.1), and a message describing the constraint. A default id is automatically generated when a new constraint is created.



4.4.4 Events

Events are used to specify discrete changes of compartment sizes, species amounts or concentrations, and parameter values. As shown below, each event is composed of an ID, an optional name, a trigger formula, an optional delay formula, and a list of event assignments. When adding a new event, a default ID is provided. The behavior of an event is that during each simulation cycle, the trigger formula is evaluated. If it was false in the previous simulation cycle, and it is now evaluating to true, the event is scheduled to occur at a time in the future specified by the delay formula or immediately if no delay formula is provided. It should be noted that since the trigger value must change from false to true, no event is scheduled if the trigger evaluates to true at the start of simulation. When an event occurs, it executes all the event assignments. Each event assignment sets a compartment size, species amount or concentration, or parameter value to the value specified by the SBML math formula (see Section 4.1) provided with the event assignment.



5 GCM Editor

The GCM editor shown below allows the user to create or modify a GCM (see Section 10). A GCM is a compact graphical representation of a genetic circuit which can later be synthesized into an SBML model. A GCM includes promoters (see Section 5.1), GCM species (see Section 5.2), influences (see Section 5.3), GCM parameters (see Section 5.4), and an optional SBML file (see Section ??). GCM species, influences, and promoters can be added, removed, or edited. Parameters can only be edited. An SBML file can also be selected to merge with the SBML generated from a GCM. This allows either customization of the SBML model or the addition of SBML constructs such as Initial Assignments, Rules, Constraints, and Events (see Section 4.4).

The screenshot shows the GCM Editor window with the title bar 'Cl_Cll.gcm'. The interface includes the following elements:

- GCM Id:** A text field containing 'Cl_Cll'.
- SBML File:** A dropdown menu showing '--none--'.
- Biochemical abstraction:** An unchecked checkbox.
- Dimerization abstraction:** An unchecked checkbox.
- List of Promoters:** An empty list box.
- List of Species:** An empty list box.
- List of Influences:** An empty list box.
- List of Parameters:** A list box containing the following parameters and their default values:
 - Activated production rate (ka), Default, .25
 - Activation binding equilibrium (Ka), Default, .0033
 - Basal production rate (kb), Default, .0001
 - Biochemical equilibrium (Kb), Default, .05
 - Degradation rate (kd), Default, .0075
 - Degree of cooperativity (nci), Default, 2
 - Dimerization equilibrium (Kd), Default, .05
 - Initial RNAP count (nr), Default, 30
 - Initial promoter count (ngi), Default, 2
- Action Buttons:**
 - Below the Promoters list: Add Promoter, Remove Promoter, Edit Promoter.
 - Below the Species list: Add Species, Remove Species, Edit Species.
 - Below the Influences list: Add Influence, Remove Influence, Edit Influence.
 - Below the Parameters list: Edit Parameter.
 - At the bottom: Save GCM, Save GCM as, Save as SBML, Save as SBML template.

To add a new item to the GCM, click on the appropriate add button. You will then be prompted to input information regarding the new item. After you have filled out the required fields, click on ok and the new item will be added into the GCM. To remove an item from the GCM, select that item and click the remove button. The item will then be removed from the GCM. However, if you try to remove species or promoters that are used in an influence, you will first have to remove the influence in order to remove the species or promoter from the model. To edit an existing item, select that item from the list and click the edit button. An editing window will open and you will be able to change the properties of that item. When you are done editing this item, click save to save the changes to the item. To merge an existing SBML file with the GCM output, click on the SBML file and select the SBML file to use. This will merge the contents of the selected SBML file with the SBML file that is generated from the GCM. Finally, there are two abstractions that can be selected to be performed when generating SBML, the biochemical and dimerization abstractions.

Once a GCM is completed, the user can save the GCM or save the GCM using a new name. The user can also save the GCM as SBML which creates an SBML file of the same name as the GCM. Finally, the user can save an SBML template which creates a blank SBML file with the same species as the GCM. This is useful for creating an SBML file which will be attached to the GCM, and includes rules, constraints, or events.

5.1 Promoters

Promoters are special species which represent the region of the DNA from which transcription is initiated. When adding or editing promoters, the user must supply a unique ID. An optional name can also be provided which is an arbitrary string description for the promoter. If desired, the user can then modify the initial promoter count (ng), the RNAP binding equilibrium (K_o), the open complex production rate (ko), the stoichiometry of production (i.e., the number of transcripts per mRNA, np), the basal production rate (kb), or the activated production rate (ka).

Dialog box titled "Promoter Editor" showing fields for editing a promoter. The "ID" field contains "PR". The "Name" field is empty. Below are several parameters with dropdown menus and input fields:

| Parameter | Value |
|------------------------------------|-------|
| Initial promoter count (ng) | 2 |
| RNAP binding equilibrium (K_o) | .033 |
| Open complex production rate (ko) | .05 |
| Stoichiometry of production (np) | 10 |
| Basal production rate (kb) | .0001 |
| Activated production rate (ka) | .25 |

Buttons: Cancel, Ok

5.2 GCM Species

GCM species are the molecules (usually proteins) produced by genes. When adding or editing a species, the user must provide a unique ID. The user can also select the type of the species to be normal, constant, or unconstrained. A normal species will result in gene production and degradation reactions being produced. A constant species will not generate any production or degradation reactions. An unconstrained species will produce a constant production and degradation reaction. The user can also specify an initial species count (ns), a Dimerization equilibrium (K_d), and a degradation rate (kd).

Dialog box titled "Species Editor" showing fields for editing a species. The "ID" field contains "CI". The "Name" field contains "The lambda repressor". The "Type" dropdown menu is set to "normal". Below are several parameters with dropdown menus and input fields:

| Parameter | Value |
|------------------------------------|-------|
| Initial species count (ns) | 0 |
| Dimerization equilibrium (K_d) | .05 |
| Degradation rate (kd) | .0075 |

Buttons: Cancel, Ok

5.3 Influences

Influences describe the relationships between the GCM species. When adding or editing an influence, the user must select an input and output species, as well as the type of influence. If the type is repression, then the input species represses the production of the output species. If the type is activation, then the input species activates the production of the output species. The user can also specify whether the influence has a promoter. If a promoter is selected, then this groups all influences using the same promoter together.

The screenshot shows the 'Influence Editor' window with the 'List of Parameters' tab selected. The parameters are as follows:

| Parameter | Value |
|-------------------------------------|-----------------------------|
| Name | CI - CII, Promoter default |
| Input | CI |
| Output | CII |
| Promoter | default |
| Type | repression |
| Biochemical | no |
| Degree of cooperativity (nc) | 2 |
| N-mer as transcription factor (nd) | 1 |
| Repression binding equilibrium (Kr) | .5 |
| Activation binding equilibrium (Ka) | .0033 |
| Biochemical equilibrium (Kb) | .05 |

For example, if there are two influences:

A -| C, Promoter P1

B -| C, Promoter P2

this will create two reactions, where in the presence of A and B, C is repressed. This would behave roughly like a NAND gate. If, on the other hand, there are two influences:

A -| C, Promoter P1

B -| C, Promoter P1

this creates one reaction, where in the presence of A or B, C is repressed. This would behave roughly like a NOR gate.

Users can also specify if the influence is a biochemical influence. A biochemical influence requires all input species belonging to the same promoter to be present in order to affect transcription.

For example, if there are two biochemical influences:

A +| C, Promoter P1

B +| C, Promoter P1


this will create two reactions, A and B combines into a complex, and the complex represses the production of C. This is a NAND gate. If, on the other hand, if biochemical is not selected, this behaves as a NOR gate.

The user can also set the value for degree of cooperativity (nc) which is the number of binding sites for transcription factors. The N-mer as transcription factor (nd) field determines how many monomers of the input species must be bound together in order to affect transcription. The user can also set the value of the repression or activation binding equilibrium (Kr and Ka). Finally, the user can specify the biochemical equilibrium (Kb).

5.4 GCM Parameters

GCM parameters are a list of global parameters that are used when generating the SBML model for the GCM. The parameter list allows the user an easy way to change all the parameter values in a convenient location. If a parameter in the GCM is set to default, it will use the value found in the GCM parameter. The GCM parameters are listed below:

| ID | Default Value | Units | Structure | Description |
|----|---------------|--------------------------------------|-----------|--------------------------------|
| nr | 30 | molecule | model | Initial RNAP count |
| ns | 0 | molecule | species | Initial species count |
| Kd | 0.05 | $\frac{1}{\text{molecule}}$ | species | Dimerization equilibrium |
| kd | 0.0075 | $\frac{1}{\text{sec}}$ | species | Degradation rate |
| ng | 2 | molecule | promoter | Initial promoter count |
| np | 10 | molecule | promoter | Stoichiometry of production |
| nc | 2 | molecule | promoter | Degree of cooperativity |
| Ko | 0.033 | $\frac{1}{\text{molecule}}$ | promoter | RNAP binding equilibrium |
| ko | 0.05 | $\frac{1}{\text{sec}}$ | promoter | Open complex production rate |
| kb | 0.0001 | $\frac{1}{\text{sec}}$ | promoter | Basal production rate |
| ka | 0.25 | $\frac{1}{\text{sec}}$ | promoter | Activated production rate |
| nd | 1 | molecule | influence | N-mer as transcription factor |
| Kr | 0.5 | $\frac{1}{\text{molecule}^{nc}}$ | influence | Repression binding equilibrium |
| Ka | 0.0033 | $\frac{1}{\text{molecule}^{(nc+1)}}$ | influence | Activation binding equilibrium |
| Kb | 0.05 | $\frac{1}{\text{molecule}}$ | influence | Biochemical equilibrium |



6 Analysis View

The analysis view is used to analyze biochemical reaction network models. The analysis view as shown below includes tabs for simulation options (see Section 6.1), abstraction options (see Section 6.2), a parameter editor (see Section 6.3), a TSD graph editor (see Section 8), and a probability graph editor (see Section 9).

The screenshot shows the 'simSBML' application window with the 'Simulation Options' tab selected. The interface includes three sets of radio buttons for configuration: 'Choose One: None (selected) | Abstraction | Logical Abstraction', 'Choose One: ODE (selected) | Monte Carlo | Markov | SBML | Network | Browser', and 'Choose One: Overwrite (selected) | Append'. Below these is a dropdown menu for 'Possible Simulators/Analyzers' with 'rkf45' selected. A 'Description Of Selected Simulator:' field shows 'Embedded Runge-Kutta-Fehlberg (4, 5) method'. Several input fields are present: 'Time Limit:' (2100.0), 'Print Interval:' (20.0), 'Maximum Time Step:' (inf), 'Absolute Error:' (1.0E-9), 'Random Seed:' (314159), 'Runs:' (1), and 'Simulation ID:' (empty). At the bottom are 'Save and Run' and 'Save Parameters' buttons.

6.1 Simulation Options

iBioSim comes with a number of simulation methods, ranging from continuous-deterministic simulation methods to discrete-stochastic simulation methods. In order to perform efficient temporal behavior analysis, various model abstraction can also be automatically applied. These routines are implemented within the reb2sac tool described in Kuwahara’s PhD Dissertation (UofUtah 2007) .

The first set of radio buttons in this tab specifies the levels of abstraction. “None” means to use no abstraction, “Abstraction” means to perform reaction-based abstraction, and “Logical Abstraction” means to perform both reaction-based and logical abstractions.

The second set of radio buttons specify the type of analysis. “ODE” is for continuous-deterministic simulation, “Monte Carlo” is for discrete-stochastic simulation, “Markov” performs temporal probability distribution analysis on finite-state Markov chain models, “sbml” outputs the model in SBML format, “Network” outputs the structure of the model in the GraphViz format for display by dot, “Browser” outputs the model in xhtml format for display in a web browser.

The last set of radio buttons asks if you want to “Overwrite” the simulation runs or if you want to “Append” more simulation runs. If you have not yet performed any simulation, this option is disabled.

The next field specifies the simulation method you want to use based on the simulation type you specified. The methods available are:

| Type | Method ID | Description |
|-------------|-----------|--|
| ODE | Euler | The forward Euler Method |
| ODE | gear1 | Gear Method M=1 |
| ODE | gear2 | Gear Method M=2 |
| ODE | rk4imp | Implicit 4th order Runge-Kutta at Gaussian points |
| ODE | rk8pd | Embedded Runge-Kutta Prince-Dormand (8,9) method |
| ODE | rkf45 | Embedded Runge-Kutta-Fehlberg (4, 5) method |
| Monte carlo | Gillespie | Gillespie's SSA direct method |
| Monte carlo | emc-sim | Monte Carlo simulation with jump count as an independent variable |
| Monte carlo | bunker | Bunker's method: next reaction time step is calculated using the average |
| Monte carlo | nmc | Monte Carlo simulation with normally distributed next reaction time |

There are some properties that need to be set for simulation. The table below specifies these:


| Field | Description |
|-------------------|--|
| Time Limit | The simulation time limit |
| Print Interval | The print time interval for each simulation run |
| Maximum Time Step | The maximum time step allowed (also minimum time step for the Euler method) |
| Absolute Error | Used by the adaptive time step ODE methods |
| Random Seed | An integer number as a seed to generate random numbers |
| Runs | The number of Monte Carlo simulation runs to perform |
| Simulation ID | Creates a simulation directory with the ID name |

6.2 Abstraction Options

This tab as shown below allows the user to set the properties of rapid equilibrium, QSSA, and operator site abstraction methods.


- Rapid Equilibrium Condition 1 specifies threshold $T1$ such that the rapid equilibrium condition fails when $T1 > E0/(S0 + k - 1/k1)$.
- Rapid Equilibrium Condition 2 specifies threshold $T2$ such that the rapid equilibrium condition fails when $T1 > k2/k - 1$.
- The QSSA condition specifies threshold T used by the QSSA abstraction method where $T > E0/(S0 + KM)$.
- The Max concentration threshold specifies the maximum number of molecules that a species can have initially and still be considered an operator site by the operator site reduction.

This tab also allows the user to select the interesting species. Interesting species are the ones that are used in the analysis, and hence are those which should never be abstracted away. This tab shows all available species, and to make a species (or set of species) interesting, highlight the species and press the Add Species button. There is also a button to remove interesting species and to clear all interesting species.



6.3 Parameter Editor

The parameter editor as shown below is similar in form to the SBML editor, but it only allows initial concentrations and parameters to be adjusted. Each of these parameters starts with the original value specified in the SBML or GCM associated with this analysis view. By changing the type to “Custom”, a new value can be entered. Changing the type back to “Original”, restores the original value. These values can also be swept by selecting the “Sweep” type. In this case, you should provide a start value, a stop value, a step amount, and a level (1 or 2). When analyzing using sweep parameters, one analysis run is produced for each value stepped through from start to stop. The parameters at level 2 are changed first. When they have all reached their stop value, the parameters at level 1 are stepped once, and the parameters at level 2 are stepped through again. This process repeats until all parameters at level 1 have stepped to their stop value.



7 Learn View

The learn view is used to discover genetic circuit connectivity from time series data. The learn view includes tabs for a data manager (see Section 7.1), a learn tool (see Section 7.2), and a TSD graph editor (see Section 8).

7.1 Data Manager

The data manager as shown below is used to both enter time series experimental data as well as bring data into the learn view. The Add button is used to create a new data file. After pressing this button, enter the name of the new data file, and then enter the data for this file using the data editor to the right. The Remove button deletes all highlighted files. Note that after highlighting one file, you can use the ctrl key to highlight additional files or the shift key to highlight a range of files. The Rename button is used to change the name of a data file. The Copy button copies a data file. The Copy From View button brings up a list of all analysis and learn views in the current project, and data from the selected view will be copied into this learn view. Finally, the Import button brings up a file browser, and it allows you to import a data file from outside this project. These files can be in time series data (TSD) format (see Section 11), comma separated value (CSV) format, or tab delimited format (DAT).

The contents of the data file highlighted on the left appear in the data editor on the right. Individual data entries can be modified, new data points can be added using the Add Data Point button, data points can be removed using the Remove Data Point button, and data points can be copied using the Copy Data Point button. When you are satisfied with all your changes, you should press the Save button to record your changes.

| | time | CI | CII |
|----------------|--------|------|------|
| sim/run-1.tsd | 0.0 | 0.0 | 0.0 |
| sim/run-10.tsd | 100.0 | 0.0 | 28.0 |
| sim/run-11.tsd | 200.0 | 74.0 | 20.0 |
| sim/run-12.tsd | 300.0 | 38.0 | 9.0 |
| sim/run-13.tsd | 400.0 | 31.0 | 4.0 |
| sim/run-14.tsd | 500.0 | 16.0 | 3.0 |
| sim/run-15.tsd | 600.0 | 32.0 | 6.0 |
| sim/run-16.tsd | 700.0 | 16.0 | 2.0 |
| sim/run-17.tsd | 800.0 | 16.0 | 23.0 |
| sim/run-18.tsd | 900.0 | 41.0 | 17.0 |
| sim/run-19.tsd | 1000.0 | 23.0 | 7.0 |
| sim/run-2.tsd | 1100.0 | 12.0 | 3.0 |
| sim/run-20.tsd | 1200.0 | 5.0 | 1.0 |
| sim/run-3.tsd | 1300.0 | 2.0 | 9.0 |
| sim/run-4.tsd | 1400.0 | 8.0 | 27.0 |
| sim/run-5.tsd | 1500.0 | 40.0 | 12.0 |
| sim/run-6.tsd | 1600.0 | 20.0 | 5.0 |
| sim/run-7.tsd | 1700.0 | 13.0 | 3.0 |
| sim/run-8.tsd | 1800.0 | 7.0 | 1.0 |
| sim/run-9.tsd | 1900.0 | 2.0 | 10.0 |
| | 2000.0 | 3.0 | 5.0 |
| | 2100.0 | 12.0 | 2.0 |

7.2 Learn Tool

The learn tool shown below uses the GeneNet algorithm described in Barker's PhD dissertation (UofUtah 2007) . To use this learn tool, adjust any options described below, if desired, then press the Save and Learn button. The resulting genetic circuit is specified using our Genetic Circuit Model (GCM) Format (see Section 10) and shown graphically using GraphViz's Dotty tool. On this tab, there are also buttons to save the parameters without learning, view the last learned circuit, save the generated circuit into the project, and view the last run log.

Below are the basic learning options as shown above are as follows:

- Minimum Number of Initial Vectors (Tn) (default=2):
Tn is a threshold value used in the CreateInfluenceVectorSet algorithm and represents the minimum number of influence vectors constructed in this algorithm.
- Maximum Influence Vector Size (Tj) (default=2):
Tj is a threshold value used in the CombineInfluenceVectors algorithm to determine the maximal size of merged influence vectors.
- Score for Empty Influence Vector (Ti) (default=0.5):
The score for an influence vector with no influences in it.
- Number of Bins (default=4):
The number of bins value specifies how many values the encoded time series data can assume.
- Equal Data Per Bins / Equal Spacing of Bins:
This radio button selects whether the auto generated levels should be determined by equaling dividing the data between the bins or by equally dividing the range of the data.
- Use Auto Generated Levels / Use User Generated Levels:
This radio button allows the user to select whether they want the levels separating the bins to be auto generated or the user would like to provide them.

- When using user provided levels, the Suggest Levels button will provide the levels that would have been auto generated as a suggestion. These levels can then edited by the user. The number of bins for each species can also be individually adjusted.

The advanced learning options shown below are as follows:

- Ratio for Activation (T_a) (default=1.15):
A probability ratio above this value results in a vote for an influence vector that has a majority of activation influences.
- Ratio for Repression (T_r) (default=0.75):
A probability ratio above this value results in a vote for an influence vector that has a majority of repression influences.
- Merge Influence Vectors Delta (T_m) (default=0.0):
Two influence vectors cannot be merged unless the difference in their scores is less than this value.
- Relax Thresholds Delta (T_t) (default=0.025):
The values of T_a and T_r are modified by this amount when these thresholds are relaxed.
- Debug Level (default=0):
This controls how much information is displayed by the GeneNet algorithm when it runs.
- Successors / Predecessors / Both (default=Successors):
This radio button selects whether successor data point pairs, predecessor data point pairs, or both are used.
- Basic FindBaseProb (default=unchecked):
When selected, the basic FindBaseProb function is used.


The screenshot shows a software window titled 'learn' with a close button. It features a tabbed interface with 'Data Manager', 'Learn' (selected), and 'TSD Graph'. Under the 'Learn' tab, there are two sub-tabs: 'Basic Options' and 'Advanced Options' (selected). The 'Advanced Options' section contains the following settings:

- Ratio For Activation (T_a): 1.15
- Ratio For Repression (T_r): 0.75
- Merge Influence Vectors Delta (T_m): 0.0
- Relax Thresholds Delta (T_t): 0.025
- Debug Level: 0 (with a spin button)
- Successors ☒ Predecessors ☐ Both ☐
- ☐ Basic FindBaseProb

At the bottom of the window, there are five buttons: 'Save and Learn' (highlighted), 'Save Parameters', 'View Circuit', 'Save Circuit', and 'View Run Log'.


8 TSD Graph Editor

The TSD graph editor appears as a tab in both analysis and learn views. TSD graphs can also be created at the top-level of the project to allow you to integrate results from several analysis or learn views. These graphs can be created using the New → TSD Graph menu option. Once created, they can be viewed and edited by clicking on the graph in the project window. An example graph is shown below.



In the TSD graph editor shown below, a graph is created by double clicking on the graph. You can then set various parameters and select what values you would like to have graphed. The parameters that you can select for a graph include:

- Title - The title of the graph.
- X-Axis Label - The label displayed for the x-axis.
- Y-Axis Label - The label displayed for the y-axis.
- X-Min - The starting value for the x-axis.
- X-Max - The ending value for the x-axis.
- X-Step - The increment for the x-axis.
- Y-Min - The starting value for the y-axis.
- Y-Max - The ending value for the y-axis.
- Y-Step - The increment for the y-Axis.
- Auto Resize Check Box - Determines whether to automatically resize the graph for best fit.



The data selection menu on the left displays all of the available sets of data that can be graphed. In particular, one can graph the average, variance, standard deviation, or results from individual simulation runs. For a top-level graph, these data sets will be organized hierarchically. Hierarchy is also introduced when simulations in an analysis view are given simulation IDs or after performing an analysis while sweeping parameter values. After selecting a data set, one can select individual species to graph and how they are to be displayed. In other words, for each species, there are the following options:

- Use Check Box - Determines whether or not this species is displayed on the graph. Checking or unchecking the box at the top changes the state for all species in the data set.
- Species Label - The name displayed in the legend.
- Drop Down Menu Of Colors - The color that is used for this species.
- Drop Down Menu Of Shapes - The shape that is used to mark the data points.
- Connect Check Box - Determines whether to connect the points with a line. Checking or unchecking the box at the top changes the state for all species in the data set.
- Visible Check Box - Determines whether shapes are visible on the line. Checking or unchecking the box at the top changes the state for all species in the data set.
- Fill Check Box - Determines whether shapes are filled on the line. Checking or unchecking the box at the top changes the state for all species in the data set.

Note that a check mark appears on a data set to indicate that some species have been selected in that data set. Also, all species can be deselected by pressing the Deselect All button.


The “Save Graph” button save the settings for the graph to a file, so when you re-open the graph, it will reload this data and display in the same way as before. The “Save As” button prompts for a filename and creates a new top-level graph with that name. Finally, the “Export” button prompts for a filename and exports the data to the given name. The extension provided for the filename is used to determine how the graph is to be exported. The supported file types are:

- csv - comma separated value data file.
- dat - column separated data file.
- eps - encapsulated postscript.
- jpg - JPEG (Joint Photographic Experts Group).
- pdf - portable document format.
- png - portable network graphics.
- svg - scalable vector graphics.
- tsd - time series data format (see Section 11).

If no extension is given, then the file type is the one specified in the file filter (default is pdf). For image (i.e., not data) file types, you will be prompted to give a desired pixel height and width for the file before the file is exported.


9 Probability Graph Editor

Probability graphs are used to display histograms for reasons that simulations terminated. This is used in conjunction with SBML constraints to determine the likelihood of various conditions. The probability graph editor appears as a tab in analysis views. Probability graphs can also be created at the top-level of the project to allow you to integrate results from several analysis views. These graphs can be created using the New → Probability Graph menu option. Once created, they can be viewed and edited by clicking on the graph in the project window. An example probability graph is shown below.



In the probability graph editor, a graph is created by double clicking on the graph. You can then set various parameters and select what values you would like to have graphed. The parameters that you can select for a graph include:

- Title - The title of the graph.
- X-Axis Label - The label displayed for the x-axis.
- Y-Axis Label - The label displayed for the y-axis.



The data selection menu on the left displays all of the available sets of data that can be graphed. For a top-level graph, these data sets will be organized hierarchically. Hierarchy is also introduced when simulations in an analysis view are given simulation IDs or after performing an analysis while sweeping parameter values. After selecting a data set, one can select individual species to graph and how they are to be displayed. In other words, for each species, there are the following options:

- Use Check Box - Determines whether or not this species is displayed on the graph. Checking or unchecking the box at the top changes the state for all species in the data set.
- Species Label - The name displayed in the legend.
- Drop Down Menu Of Colors - The color that is used for this species.

Note that a check mark appears on a data set to indicate that some species have been selected in that data set. Also, all species can be deselected by pressing the Deselect All button.

The “Save Graph” button save the settings for the graph to a file, so when you re-open the graph, it will reload this data and display in the same way as before. The “Save As” button prompts for a filename and creates a new top-level graph with that name. Finally, the “Export” button prompts for a filename and exports the data to the given name. The extension provided for the filename is used to determine how the graph is to be exported. The supported file types are:

- eps - encapsulated postscript.
- jpg - JPEG (Joint Photographic Experts Group).
- pdf - portable document format.
- png - portable network graphics.

- svg - scalable vector graphics.

If no extension is given, then the file type is the one specified in the file filter (default is pdf). For image (i.e., not data) file types, you will be prompted to give a desired pixel height and width for the file before the file is exported.

10 Genetic Circuit Model Format

Our genetic circuit model (gcm) format specifies a genetic circuit using the same format used by the GraphViz graph drawing tool. The vertices in the graph are the species in the genetic circuit, and the edges in the graph represent the activation and repression relationships between the species. An activation relationship is shown with a blue (blue4) arrow (vee) and a repression relationship is shown with a red (firebrick4) tee. The label field in the species declaration is the name of the species. The arrowhead field in the relationship declaration represents the type of relationship between the species. Repression is labeled with a tee and activation is labeled with a vee. The label field in the relationship declaration represents how many molecules are necessary to activate or repress the production of the species. An example is shown below for a simple genetic circuit in which the species CI represses CII while CII activates CI production. The $s1 \rightarrow s2$ edge has a label field of “2” which means two molecules of CI are required to form a dimer to repress CII.

```
digraph G {
  s1 [shape=ellipse,color=black,label="CI"];
  s2 [shape=ellipse,color=black,label="CII"];
  s2 -> s1 [color="blue4",arrowhead=vee];
  s1 -> s2 [color="firebrick4",label="2",arrowhead=tee];
}
```

More advanced behavior can be modeled by using extra fields. The promoter field groups a set of species together. The examples below shows how the promoter field works. In the genetic circuit model below, species A represses the production of species B and C, independently. If there was exactly 1 molecule of species A, it would only be able to repress production of species B or C, but not both.

```
digraph G {
  s1 [shape=ellipse,color=black,label="A"];
  s2 [shape=ellipse,color=black,label="B"];
  s3 [shape=ellipse,color=black,label="C"];
  s1 -> s2 [color="blue4",arrowhead=tee];
  s1 -> s3 [color="blue4",arrowhead=tee];
}
```

With the promoter field, one species A now represses the promoter “P1”, which produces both species B and C. This means that one molecule of species A will repress the production of both species B and C.

```
digraph G {
  s1 [shape=ellipse,color=black,label="A"];
  s2 [shape=ellipse,color=black,label="B"];
  s3 [shape=ellipse,color=black,label="C"];
```

```

s1 -> s2 [color="blue4",arrowhead=tee,promoter="P1"];
s1 -> s3 [color="blue4",arrowhead=tee,promoter="P1"];
}

```

The promoter field can also be used to separate production reactions. In the example below, both species A and B can repress the production of species C. If either is present, then very little C will be produced. This behavior is like a NOR gate.

```

digraph G {
  s1 [shape=ellipse,color=black,label="A"];
  s2 [shape=ellipse,color=black,label="B"];
  s3 [shape=ellipse,color=black,label="C"];
  s1 -> s3 [color="blue4",arrowhead=tee];
  s2 -> s3 [color="blue4",arrowhead=tee];
}

```

However, if there needs to be two different sources of production for species C, the promoter field can be used to accomplish this. In the example below, A represses the production of C by binding to the P1 promoter, and B represses the production of C by binding to the P2 promoter. Both A and B need to be present to fully repress the level of C. If either is at a low level, then the level of C will be high. This behavior is like a NAND gate.

```

digraph G {
  s1 [shape=ellipse,color=black,label="A"];
  s2 [shape=ellipse,color=black,label="B"];
  s3 [shape=ellipse,color=black,label="C"];
  s1 -> s3 [color="blue4",arrowhead=tee,promoter="P1"];
  s2 -> s3 [color="blue4",arrowhead=tee,promoter="P2"];
}

```

The example below shows how to model an AND gate. The species A and B have the constant flag set to true. This means that A and B have no production and degradation reactions. The reactions contain a promoter label “P1”. This means that the species C can be activated by both s1 and s2. Combined with the type flag of biochemical, this creates a biochemical reaction where species A and B combine together to form a complex to activate production of species C.

```

digraph G {
  s1 [shape=ellipse,color=black,label="A",const=true];
  s2 [shape=ellipse,color=black,label="B",const=true];
  s3 [shape=ellipse,color=black,label="C"];
  s1 -> s3 [color="blue4",arrowhead=vee,promoter="P1",type=biochemical];
  s2 -> s3 [color="blue4",arrowhead=vee,promoter="P1",type=biochemical];
}

```

11 Time Series Data Format

The time series data (tsd) format is composed of a parenthesized and comma-separated set of time points. Each time point is composed of a parenthesized and comma-separated set of data for that

time point. This first time point is composed of a set of strings that are the labels for the data entries. The first entry in each time point is by convention the time for that time point. Below is an example simulation of the species *CI* and *CII* from 0 to 1000 seconds with time points separated by 100 seconds.


(("time","CI","CII"), (0,0,0), (100,0,19), (200,20,25), (300,19,18), (400,17,20), (500,17,46), (600,26,40), (700,43,43), (800,63,28), (900,72,34), (1000,72,28))

12 Tutorial

The example described in this section constructs a simple model for the *cI* and *cII* genes and the P_R and P_{RE} promoters from the phage λ decision circuit. It illustrates many of the features of iBioSim. After starting iBioSim, complete the following steps to create an SBML model for this example:

1. Select **File** \rightarrow **New** \rightarrow **Project**
2. Browse to desired path and create a project named `lambda`.
3. Select **File** \rightarrow **New** \rightarrow **SBML Model**
4. Enter `lambda` which will open the SBML editor.
5. Select **Definitions/Types** tab.
6. Select **Add Unit** and enter `per_second` as the ID.
7. Select **Add to List**, select `second` as the kind, change the exponent to -1 , and click **Add**.
8. Click **Add** in the **Unit Definition Editor**.
9. Select **Add Unit** and enter `per_second_mole` as the ID.
10. Select **Add to List**, select `second` as the kind, change the exponent to -1 , and click **Add**.
11. Select **Add to List**, select `mole` as the kind, change the exponent to -1 , and click **Add**.
12. Click **Add** in the **Unit Definition Editor**.
13. Select **Main Elements** tab.
14. Select `default` in the list of Compartments, and press the **Edit Compartment** button, and change the compartment ID to `cell` and size to 1.
15. Select **Add Species** and enter `CI` as the ID.
16. Select **Add Species** and enter `CI2` as the ID.
17. Select **Add Reaction** and enter `Dimerize` as the ID and change reversible to `true`.
18. Select **Add Reactant** and select `CI` as the species and change the stoichiometry to 2.
19. Select **Add Product** and select `CI2` as the species.
20. Highlight `kf` and select **Edit Selected Parameter**, change `kf` to `k2f`, and change the units to `per_second_mole`.

21. Highlight `kr` and select **Edit Selected Parameter**, change `kr` to `k2r`, change the value to 1.0, and change the units to `per_second`.
22. Select **Use Mass Action**.
23. Select **Add**.
24. Select **Save**.
25. Highlight `lambda.sbml`, using right mouse button, select **View Network**.
26. Highlight `lambda.sbml`, using right mouse button, select **View in Browser**.
27. Go back to the SBML editor complete the construction of the chemical reaction network shown below:



| Constant | Value | Constant | Value | Constant | Value | Constant | Value |
|--------------|---------------|--------------|--------------------------|--------------|---------------------------|-------------|-------------------------|
| <i>KPRE2</i> | 0.01 M^{-1} | <i>KPRE4</i> | 0.00161 M^{-2} | <i>kPREb</i> | 0.00004 sec^{-1} | <i>kPRE</i> | 0.015 sec^{-1} |
| <i>n</i> | 10 | <i>KOR9</i> | 0.69422 M^{-1} | <i>KOR10</i> | 0.06568 M^{-2} | <i>kPR</i> | 0.014 sec^{-1} |
| <i>K2</i> | 0.1 M^{-1} | <i>k1</i> | 0.0007 sec^{-1} | <i>k10</i> | 0.002 sec^{-1} | | |

Set an initial amount of 1.0 for PRE and OR, 30.0 for RNAP, and 0.0 for the rest.

28. Highlight `lambda.sbml`, using the right mouse button, select **Create Analysis View**, and enter the name `sim`.
29. In the newly opened window, select **Monte Carlo**.
30. Also, in this window, change the time limit to 2100.0, print interval to 100.0, and runs to 20.
31. Finally, select **Save and Run** at the bottom of the window.
32. After the simulation completes (it may take a little while), click on the graph tab.
33. Click on the graph to bring up the graph editor. Highlight Average, if not already highlighted, select CI2 and CII, change Title to “Average”, change X-Axis Label to “Time (seconds)”, and change Y-Axis Label to “Number of Molecules”. Press the OK button. Click on Export and enter file name of `average.jpg`.

34. Click on the graph to bring up the graph editor. Press the “Deselect All” button. Highlight run-1, select CI2 and CII, change Title to “Run-1”, change X-Axis Label to “Time (seconds)”, and change Y-Axis Label to “Number of Molecules”. Press the OK button. Click on Export and enter file name of `run1.jpg`.
35. Click on the graph to bring up the graph editor. Press the “Deselect All” button. Highlight Standard Deviation, select CI2 and CII, change Title to “Standard Deviation”, change X-Axis Label to “Time (seconds)”, and change Y-Axis Label to “Number of Molecules”. Press the OK button. Click on Export and enter file name of `stddev.jpg`.
36. Click on the SBML editor tab. Change the initial amounts of OR and PRE to 10. Press the Save and Run button.
37. Click on the graph tab and following the steps above, create the following plots `average_10.jpg`, `run1_10.jpg`, and `stddev_10.jpg`. How are these different than the first set of plots?
38. Simulate your lambda model with BioSim using the ODE method rkf45 with a time limit of 2100 and print interval of 50. Make a note of the simulation time and plot CI2 and CII. Next, simulate using the Euler method. Make a note of the simulation time and add CI2 and CII from the euler results to your graph. How do the simulation times and results compare? Change the time step and rerun the Euler method. Repeat until the results match up well. What time step is required for a good match? How do the simulation times compare?
39. Open the lambda model in an SBML editor, and add a new species “CI_t” with a 0 initial amount.
40. Click on the “Initial Assignments/Rules/Constraints/Events” tab and press the “Add Rule” button.
41. Select Assignment Type, select Variable “CI_t”, and enter the following as the Rule:

$$2 * CI2 + CI$$
42. Press the Save SBML button.
43. Create or open an analysis view on your lambda model.
44. Select “None” and “Network”, and press the Save and Run button. Count the number of species and reactions in your model.
45. Select “Monte Carlo”, set the Time Limit to 2100.0, Print Interval to 100.0, and Runs to 20. Press the Save and Run button and record the simulation time.
46. Create a new analysis view for your lambda model and select “Abstraction” this time.
47. In the abstraction tab, change rapid equilibrium conditions 1 and 2 to 1000.0 as well as QSSA condition to 1000.0. Select “Network” then save and run. Count the number of species and reactions in your model.
48. Select “Monte Carlo”, set the Time Limit to 2100.0, Print Interval to 100.0, and Runs to 20. Press the Save and Run button and record the simulation time.
49. Select File → New → Graph and enter a name for your new top-level graph.

50. Click on the graph and find the average simulation results from your original model and graph CII and CIIt. Also, add to this graph from your abstracted model CII and CI. How well do they compare? Send me by email a jpg of your result.
51. Go back to your analysis view in which you did abstraction and change all the conditions back to 0.1. Regenerate the network and record the number of species and reactions. Regenerate the simulation results and record the simulation time.
52. Go back to the top-level graph which should have updated results. How does it compare now?
53. Create a new project named “xor”, select File → New → Genetic Circuit Model, and enter the name **xor**.
54. Add species *A* and *B* of type constant. Be sure to fill in both the ID and Name fields.
55. Add species *Abar*, *Bbar*, *X*, *Y*, and *C* of type normal. Again, fill in the ID and Name fields.
56. Add an influence of type repression with input *A* and output *Abar* and dimer of 2. You have just created an inverter. Create another inverter from *B* to *Bbar*.
57. Add a promoter named PX1 and another promoter named PX2.
58. Add a repression influence from *A* to *X* and use the promoter PX1 with dimer of 2.
59. Add a repression influence from *Bbar* to *X* and use the promoter PX2 with dimer of 2. You have just created a NAND gate.
60. Create another NAND gate with inputs *Abar* and *B* and output *Y*.
61. Finally, create another NAND gate with inputs *X* and *Y* and output *C*.
62. Change the **decay** parameter to 0.01.
63. Press the “Save GCM” and “Save as SBML” buttons.
64. Create an analysis view for your **xor.sbml** file.
65. Select the User Defined Data tab and click on the “Use User Defined Data” button.
66. Add a data point for A at time step 1000 to go to 20.
67. Add a data point for B at time step 2000 to go to 20.
68. Add a data point for A at time step 3000 to go to 0.
69. Add a data point for B at time step 4000 to go to 0.
70. Select the options tab, select “abstraction”, and set a time limit of 5000, a print interval of 200, and 20 runs.
71. Select the abstraction tab, and add species *C* to the interesting species list.
72. Press the “Save Parameters” button then “Save and Run” button.
73. Create a graph that includes *A*, *B*, and *C*. Email a jpg of this graph to me.

74. Close your analysis view and edit your `xor.gcm`. Try adjusting some of the global parameters. Remember to save both your `gcm sbml` files. Reopen your analysis view and re-run your simulation. Send me a few different graphs. Change the titles of the graphs and provide a description in your email of the graphs as to what you changed.

13 Reporting Bugs and Feature Requests

In order to report a bug or to request a change or feature, please send an email to:

`atacs-bugs@vlsigroup.ece.utah.edu`.

The subject line must begin with one of the following keywords or the mail will be filtered by our spam filters:

- BUG - error or crash of the software
- CHANGE - something which can be improved
- FEATURE - something new

14 Credits

The iBioSim tool is being developed at the University of Utah by Chris Myers , Nathan Barker , Hiroyuki Kuwahara , Curtis Madsen , and Nam Nguyen Nathan Barker is now with Southern Utah University, and Hiroyuki Kuwahara is now with the Centre for Computational and System Biology in Trento, Italy.